



EDITION 2018

Iris



BIOTECH GMBH



BIOTINYLATION

Empowering Peptide Innovation

Headquarters:

Iris Biotech GmbH
Adalbert-Zoellner-Str. 1
95615 Marktredwitz, Germany
Phone: +49 9231 97121 - 0
Fax: +49 9231 97121 - 99
Email: info@iris-biotech.de
Internet: www.iris-biotech.de

Office Belgium:

Phone: +49 9231 97121 - 28
Fax: +49 9231 97121 - 99

Distribution Partners:

India:

Sumit Biosciences Pvt Ltd.
D Wing, Krishna Complex
Subhash Road - A, Vile Parle (East)
Mumbai - 400 057, India
Phone: +91 22 6693 8885 / 6
Fax: +91 22 6699 0665
Skype: [sumit.bio](https://www.skype.com/name/sumit.bio)
Email: sumit_exports@yahoo.com
Internet: www.sumitbiomedical.com
www.sumitglobal.com

Japan:

BizCom Japan, Inc.
7F K2 buiding 5-28-1, Higashigotanda, Shinagawa-ku,
Tokyo 141-0022, Japan
Phone: +81 3 6277 3233
Fax: +81 3 6277 3265
Mobile: +81 70 5546 6277
Email: irisbiotech@bizcomjapan.co.jp

USA & Canada:

Peptide Solution LLC
8110 S. Houghton Road, Suite 158 PMB 152
Tucson, AZ 85747, USA
Phone: +1 502 314 1515
Fax: +1 520 422 3333
Email: sales@easypeptide.com
Internet: www.easypeptide.com

Content


1. Biotinylation Reagents	4
Reagents for Tyrosine- / Protein Biotinylation.....	4
Reversible Biotinylation / Reversible Biotin Binding	5
(Classical) Biotinylation Reagents	5
Biotin Click Reagents	6
Biotin with Amino Acids	7
PEG Containing Biotinylation Reagents.....	9
Biotin-PEG-Acids / Amine Reactive Biotinylation Reagents.....	9
Biotin-PEG-Amines	12
Biotin-PEG-Maleimides	14
Biotin-PEG-Thiols.....	15
Other Biotin-PEG Reagents	16
2. Terms and Conditions of Sales.....	17





1. Biotinylation Reagents

The affinity and specificity of the avidin-biotin interaction have been exploited for numerous applications in immunology, histochemistry, and affinity chromatography, to name a few. Biotinylation is a common technique for transforming poorly detectable molecules into probes that can be recognized by labeled biotin-binding proteins or an affinity capture matrix. Antibodies can be "tagged" with biotinylation reagents and used to probe cells or tissues for specific antigens or haptens.

Reagents for Tyrosine- / Protein-Biotinylation

For long time it is known that tyramine compounds are converted to highly reactive radicals by horseradish peroxidase in presence of H_2O_2 . In vicinity of proteins, these radicals preferentially react with surface exposed tyrosines. Therefore, biotin tyramide is a perfect reagent for protein biotinylation. Numerous applications have been developed and commercialized using this excellent reactivity of biotin tyramide for signal amplification and sensitive detection of very small concentrations. Biotin tyramide and also a variation with a short PEG spacer in between are available to improve flexibility and solubility.



		Article No.	Quantity	Price
LS-3500	Biotin Tyramide			
		(3aS,4S,6aR)-hexahydro-N-[2-(4-hydroxyphenyl)ethyl]-2-oxo-1H-thieno[3,4-d]imidazole-4-pentanamide CAS-NO: 41994-02-09 FORMULA: $C_{18}H_{22}N_4O_3S$ MOL WEIGHT: 363,47 g/mol		LS-3500.0250 250 mg € 150,00 LS-3500.1000 1 g € 375,00 LS-3500.5000 5 g € 1.500,00
LS-3490	Biotin-AEEA-Phenol			
		N-(2-(2-(2-(4-hydroxyphenethylamino)-2-oxoethoxy)ethoxy)ethyl)-5-((3aS,4S,6aR)-2-oxo-hexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide FORMULA: $C_{24}H_{36}N_4O_6S$ MOL WEIGHT: 508,63 g/mol		LS-3490.0100 100 mg € 275,00 LS-3490.0250 250 mg € 550,00 LS-3490.1000 1 g € 1.750,00
LS-3570	Biotin-SS-Tyramide			
		N-(2-(3-(4-hydroxyphenethylamino)-3-oxopropyl)disulfanyl)ethyl)-5-(2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)pentanamide CAS-NO: 678975-20-7 FORMULA: $C_{23}H_{34}N_4O_4S_3$ MOL WEIGHT: 526,74 g/mol		LS-3570.0250 250 mg € 350,00 LS-3570.0001 1 g € 1.000,00 LS-3570.0005 5 g € 4.000,00
Reversible biotinylation by cleavable linker.				
LS-1660	Desthiobiotin-Tyramide			
		N-(4-hydroxyphenethyl)-6-((4R,5S)-5-methyl-2-oxoimidazolidin-4-yl)hexanamide FORMULA: $C_{18}H_{22}N_4O_3$ MOL WEIGHT: 333,43 g/mol SYNONYMS: Desthiobiotin Phenol		RL-2950.0000 please inquire
Biotin-Tyramide analog that binds less tightly to biotin-binding proteins. It can be displaced competitively by biotin. Streptavidin-based ligands can be gently stripped from desthiobiotin-tyramide-labeled targets with buffered biotin solutions.				

References:


- > The Oxidation of Tyramine, Tyrosine, and Related Compounds by Peroxidase; A.J. Gross, I.W. Sizer; *J. Biol. Chem.* 1959; **234**: 1611-1614.
- > Catalyzed reporter deposition, a novel method of signal amplification application to immunoassays; M.N. Bobrow et al.; *Journal of Immunological Methods* 1989; **125**: 279-285. doi: 10.1016/0022-1759(89)90104-X.
- > Catalyzed reporter deposition, a novel method of signal amplification: II. Application to membrane immunoassays; M.N. Bobrow et al.; *Journal of Immunological Methods* 1991; **137**: 103-112. doi: 10.1016/0022-1759(91)90399-Z.
- > Tyramide signal amplification for analysis of kinase activity by intracellular flow cytometry; M.R. Clutter et al.; *Cytometry A*. 2010; **77** (11): 1020-31; doi: 10.1002/cyto.a.20970.
- > Proteomic mapping of mitochondria in living cells via spatially restricted enzymatic tagging; H.W. Rhee et al., *Science* 2013; **339**: 1328-31. doi: 10.1126/science.1230593.
- > WO2008128352 A1.

Reversible Biotinylation / Reversible Biotin Binding

In some cases, reversible biotin-mediated binding is desired. Therefore, we have biotinylating variants with cleavable linkers or less strongly binding biotin analogues.

The usually almost irreversible binding between biotin and streptavidin can be weakened by using Desthiobiotin instead of Biotin. Examples in our portfolio are D-Desthiobiotin (LS-1650; below) and Desthiobiotin-Tyramide (LS-1660; left page). Other derivatives are feasible on custom synthesis basis. Desthiobiotin can be displaced competitively from streptavidin by biotin. Streptavidin-based ligands can be gently stripped from desthiobiotin-labeled targets with buffered biotin solutions avoiding harsh denaturation conditions such as 6 M GdmCl.





Alternatively, Biotin-SS-Tyramide (LS-3570; left page) and Biotin-dPEG(4)-S-S-NHS (PEG1910; below) are disulfide bridge containing biotinylating constructs that can be cleaved using reducing agents such as glutathione or by reducing conditions inside certain cell compartments. One possible application of Biotin-SS-Tyramide is the targeted biotinylation of proteins and protein clusters on cell surfaces using HRP-labeled antibodies. Labeled proteins are then purified by using streptavidin and released by reductive cleavage of the linker disulfide bond.







References:

- > Easily reversible desthiobiotin binding to streptavidin, avidin, and other biotin-binding proteins: uses for protein labeling, detection, and isolation; J.D. Hirsch et al.; *Analytical Biochemistry* 2002; **308** (2): 343-357. doi:10.1016/S0003-2697(02)00201-4
- > A novel tandem affinity purification strategy for the efficient isolation and characterization of native protein complexes; C.J. Gloeckner et al.; *Proteomics* 2007; **7** (23): 4228-4234. doi:10.1002/pmic.200700038
- > New Insights into the DT40 B Cell Receptor Cluster Using a Proteomic Proximity Labeling Assay; X.-W. Li, J. S. Rees, P. Xue, H. Zhang, S. W. Hamaia, B. Sanderson, P. E. Funk, R. W. Farndale, K. S. Lilley, S. Perrett and A. P. Jackson; *Journal of Biological Chemistry* 2014; **289**: 14434-14447. doi:10.1074/jbc.M113.529578
- > A Catalysis-Based Selection for Peroxidase Antibodies with Increased Activity; J. Yin, J. H. Mills and P. G. Schultz; *Journal of the American Chemical Society* 2004; **126**: 3006-3007. doi:10.1021/ja039198o






(Classical) Biotinylating Reagents

			Article No.	Quantity	Price
LS-1070 D-Biotin Hexahydro-2-oxo-1H-thieno[3,4-d]imidazole-4-pentanoic acid CAS-NO: 58-85-5 FORMULA: C ₁₀ H ₁₆ N ₂ O ₃ MOL WEIGHT: 244,31 g/mol		LS-1070.0005	5 g	€ 35,00	
		LS-1070.0025	25 g	€ 125,00	
		LS-1070.0100	100 g	€ 400,00	
RL-1006 Biotin-OSu N-Hydroxysuccinimidobiotin CAS-NO: 35013-72-0 FORMULA: C ₁₄ H ₁₉ N ₃ O ₅ MOL WEIGHT: 341,39 g/mol		RL-1006.0001	1 g	€ 85,00	
		RL-1006.0005	5 g	€ 300,00	
		RL-1006.0025	25 g	€ 1.200,00	
PEG1910 Biotin-dPEG(4)-S-S-NHS 1-Biotinamino-15-oxo-3,6,9,12-tetraoxa-19,20-dithia-16-azatricosan-23-oic acid succinimidyl ester CAS-NO: 1260247-51-5 FORMULA: C ₃₀ H ₄₉ N ₃ O ₁₁ S ₃ MOL WEIGHT: 751,93 g/mol		PEG1910.0050	50 mg	€ 300,00	
		PEG1910.0500	500 mg	€ 1.150,00	
Reversible biotinylation by cleavable linker.					
LS-1650 D-Desthiobiotin (4R,5S)-5-methyl-2-oxo-4-imidazolidinehexanoic acid CAS-NO: 533-48-2 FORMULA: C ₁₀ H ₁₈ N ₂ O ₃ MOL WEIGHT: 214,27 g/mol SYNONYMS: (+)-Desthiobiotin; 5-Methyl-2-oxo-4-imidazolinecaproic acid		RL-2950.0000	please inquire		

Biotin analog that binds less tightly to biotin-binding proteins. It can be displaced competitively by biotin. Streptavidin-based ligands can be gently stripped from desthiobiotin-labeled targets with buffered biotin solutions.

		Article No.	Quantity	Price
RL-8615 Biotin-PFP Biotin pentafluorophenyl ester CAS-NO: 120550-35-8 FORMULA: C ₁₆ H ₁₅ F ₅ N ₂ O ₃ S MOL WEIGHT: 410,36 g/mol		RL-8615.0001	1 g	€ 95,00
		RL-8615.0005	5 g	€ 350,00
		RL-8615.0025	25 g	€ 1.400,00
		Explanation see p.9 (Amine Reactive Biotinylating Reagents)		
RL-2025 Biotin-Ahx-OH 6-Biotinylamino-hexanoic acid CAS-NO: 72040-64-3 FORMULA: C ₁₆ H ₂₇ N ₃ O ₄ S MOL WEIGHT: 357,47 g/mol		RL-2025.0001	1 g	€ 110,00
		RL-2025.0005	5 g	€ 400,00
		RL-2025.0025	25 g	€ 1.600,00
		Explanation see p.9 (Amine Reactive Biotinylating Reagents)		
RL-2020 Biotin-Ahx-NHS 6-Biotinylamino-hexanoic acid-N-hydroxysuccinimidyl ester CAS-NO: 72040-63-2 FORMULA: C ₂₀ H ₃₀ N ₄ O ₅ S MOL WEIGHT: 454,54 g/mol		RL-2020.0500	500 mg	€ 126,00
		RL-2020.0001	1 g	€ 196,00
		RL-2020.0005	5 g	€ 700,00
		RL-2020.0025	25 g	€ 2.800,00
RL-2030 Biotin cadaverine nTFA 5-(Biotinamido)pentylamine nTFA CAS-NO: 115416-38-1 net FORMULA: C ₁₅ H ₂₈ N ₄ O ₂ S nC ₂ HF ₃ O ₂ MOL WEIGHT: 328,47 net g/mol		RL-2030.0100	100 mg	€ 84,00
		RL-2030.0250	250 mg	€ 140,00
		RL-2030.0500	500 mg	€ 252,00
		RL-2030.0001	1 g	€ 392,00
		RL-2030.0005	5 g	€ 1.400,00
		Explanation see p.9 (Amine Reactive Biotinylating Reagents)		







Biotin Click Reagents

		Article No.	Quantity	Price
PEG4950 Biotin-PEG(4)-alkyne 15-[D(+)-Biotinylamino]-4,7,10,13-tetraoxapentadec-1-yne CAS-NO: 1262681-31-1 FORMULA: C ₂₁ H ₃₈ N ₃ O ₅ S MOL WEIGHT: 457,58 g/mol		PEG4950.0250	250 mg	€ 350,00
		PEG4950.0001	1 g	€ 1.150,00
		Explanation see p.9 (Amine Reactive Biotinylating Reagents)		
		Explanation see p.9 (Amine Reactive Biotinylating Reagents)		
PEG4940 Biotin-PEG(3)-N₃ 11-[D(+)-Biotinylamino]-1-azido-3,6,9-trioxaundecane CAS-NO: 875770-34-6 FORMULA: C ₁₈ H ₃₂ N ₃ O ₅ S MOL WEIGHT: 444,55 g/mol		PEG4940.0000	please inquire	
		Explanation see p.9 (Amine Reactive Biotinylating Reagents)		
PEG7800 Biotin-PEG(3)-ATFB alpha-Biotin-omega-(4-azido-2,3,5,6-tetrafluorobenzoylamidyl)-tri(ethylene glycol) FORMULA: C ₂₇ H ₃₇ F ₄ N ₄ O ₆ S MOL WEIGHT: 663,68 g/mol		PEG7800.0025	25 mg	€ 250,00
		PEG7800.0100	100 mg	€ 425,00
		Aryl azide that can be used for both click chemistry and for photoactivated labeling. For further information see www.iris-biotech.de/peg7800		
PEG2065 Biotin-TEG-ATFB Biotin-triethylenglycol-(p-azido-tetrafluorobenzamide) CAS-NO: 1264662-85-2 FORMULA: C ₂₇ H ₃₇ F ₄ N ₄ O ₅ S MOL WEIGHT: 663,68 g/mol		PEG2065.0025	25 mg	€ 250,00
		PEG2065.0100	100 mg	€ 425,00
		Aryl azide that can be used for both click chemistry and for photoactivated labeling. For further information see www.iris-biotech.de/peg2065		
PEG4330 Biotin-dPEG™(7)-N₃ alpha-Biotin-omega-azido hepta(ethylene glycol) CAS-NO: 1334172-75-6 FORMULA: C ₂₆ H ₄₈ N ₃ O ₅ S MOL WEIGHT: 620,76 g/mol		PEG4330.0100	100 mg	€ 325,00
		PEG4330.1000	1 g	€ 1.750,00
		Explanation see p.9 (Amine Reactive Biotinylating Reagents)		

		Article No.	Quantity	Price
PEG4340	Biotin-dPEG™(11)-N₃ [2-(2-aminoethoxy)ethoxy]acetic acid tert-butyl ester*HCl CAS-NO: 956494-20-5 FORMULA: C ₃₄ H ₆₄ N ₆ O ₁₃ S MOL WEIGHT: 796,97 g/mol	PEG4340.0100	100 mg	€ 380,00
		PEG4340.1000	1 g	€ 1.800,00
PEG4350	Biotin-dPEG™(23)-N₃ alpha-Biotin-omega-azido 23(ethylene glycol) CAS-NO: 956494-20-5 FORMULA: C ₅₈ H ₁₁₀ N ₆ O ₂₅ S MOL WEIGHT: 1325,6 g/mol	PEG4350.0100	100 mg	€ 440,00
		PEG4350.1000	1 g	€ 2.050,00
RL-2520	DBCO-PEG(4)-Biotin Dibenzoazacyclooctyne-tetra(ethylene glycol)-biotin CAS-NO: 125942-07-4 FORMULA: C ₃₉ H ₅₁ N ₅ O ₈ S MOL WEIGHT: 749,92 g/mol	RL-2520.0010	10 mg	€ 230,00
		RL-2520.0025	25 mg	€ 390,00
		RL-2520.0100	100 mg	€ 1.120,00
PEG6820	Biotin-PEG(12)-DBCO Biotinyl-dodeca(ethylene glycol)-amido-dibenzoazacyclooctyne FORMULA: C ₅₅ H ₈₃ N ₅ O ₁₆ S MOL WEIGHT: 1102,34 g/mol	PEG6820.0010	10 mg	€ 325,00
		PEG6820.0025	25 mg	€ 500,00
		PEG6820.0100	100 mg	€ 1.600,00

Biotin with Amino Acids

		Article No.	Quantity	Price
FAA3720	Fmoc-L-Cys(Biotin)-OH N-alpha-(9-Fluorenylmethoxycarbonyl)-S-biotinyl-L-cysteine FORMULA: C ₂₈ H ₃₁ N ₃ O ₆ S ₂ MOL WEIGHT: 569,69 g/mol	FAA3720.0250	250 mg	€ 90,00
		FAA3720.0500	500 mg	€ 162,00
		FAA3720.0001	1 g	€ 252,00
		FAA3720.0005	5 g	€ 900,00
		FAA3720.0025	25 g	€ 3.600,00
LS-3510	Biocytin N-epsilon-Biotinyl-L-Lysine CAS-NO: 576-19-2 FORMULA: C ₁₆ H ₂₈ N ₄ O ₄ S MOL WEIGHT: 372,48 g/mol	LS-3510.0250	250 mg	€ 95,00
		LS-3510.1000	1 g	€ 375,00
		LS-3510.5000	5 g	€ 1.500,00
Biocytin is a multi-purpose marker for neuroanatomical studies and the ideal substrate for biotinidase assays. Its strong complex formation with avidin and streptavidin conjugates opens up a broad range of detection methods.				
HAA3430	H-L-Lys(Biotin)-NH₂ N-epsilon-biotin-L-lysine amide CAS-NO: 61125-53-9 FORMULA: C ₁₆ H ₂₈ N ₄ O ₃ S MOL WEIGHT: 371,5 g/mol	HAA3430.0100	100 mg	€ 170,00
		HAA3430.0250	250 mg	€ 280,00
		HAA3430.0500	500 mg	€ 510,00
		HAA3430.1000	1 g	€ 785,00
		HAA3430.5000	5 g	€ 2.800,00
BAA1276	Boc-L-Lys(Biotin)-OH N-alpha-t-Butyloxycarbonyl-N-epsilon-biotinyl-L-lysine CAS-NO: 62062-43-5 FORMULA: C ₂₁ H ₃₆ N ₄ O ₆ S MOL WEIGHT: 472,6 g/mol	BAA1276.0500	500 mg	€ 81,00
		BAA1276.0001	1 g	€ 126,00
		BAA1276.0005	5 g	€ 450,00
		BAA1276.0025	25 g	€ 1.800,00
BAA1038	Boc-D-Lys(Biotin)-OH N-alpha-t-Butyloxycarbonyl-N-epsilon-(Biotin)-D-lysine FORMULA: C ₂₁ H ₃₆ N ₄ O ₆ S MOL WEIGHT: 472,61 g/mol	BAA1038.0500	500 mg	€ 108,00
		BAA1038.0001	1 g	€ 168,00
		BAA1038.0005	5 g	€ 600,00
		BAA1038.0025	25 g	€ 2.400,00

		Article No.	Quantity	Price
FAA1443 Fmoc-L-Lys(Biotin)-OH N-alpha-(9-Fluorenylmethoxycarbonyl)-N-epsilon-biotinyl-L-lysine CAS-NO: 146987-10-2 FORMULA: C ₃₁ H ₃₈ N ₄ O ₆ S MOL WEIGHT: 594,7 g/mol		FAA1443.0500	500 mg	€ 90,00
		FAA1443.0001	1 g	€ 140,00
		FAA1443.0005	5 g	€ 500,00
		FAA1443.0025	25 g	€ 2.000,00
FAA1451 Fmoc-D-Lys(Biotin)-OH N-alpha-(9-Fluorenylmethoxycarbonyl)-N-epsilon-biotinyl-D-lysine CAS-NO: 110990-09-5 FORMULA: C ₃₁ H ₃₈ N ₄ O ₆ S MOL WEIGHT: 594,7		FAA1451.0250	250 mg	€ 100,00
		FAA1451.0500	500 mg	€ 180,00
		FAA1451.0001	1 g	€ 280,00
		FAA1451.0005	5 g	€ 1.000,00
FAA4670 Fmoc-L-Lys(Biotin-Ahx)-OH N-alpha-(9-Fluorenylmethoxycarbonyl)-N-epsilon-[6-(biotinylamino)hexanoyl]-L-lysine CAS-NO: 160158-05-4 FORMULA: C ₃₇ H ₄₉ N ₅ O ₇ S MOL WEIGHT: 707,88 g/mol		FAA4670.0250	250 mg	€ 120,00
		FAA4670.0500	500 mg	€ 220,00
		FAA4670.0001	1 g	€ 340,00
		FAA4670.0005	5 g	€ 1.200,00
PEG4440 Fmoc-L-Lys(dPEG™(4)-Biotin)-OH N-alpha-(9-Fluorenylmethoxycarbonyl)-N-epsilon-[15-(biotinamido)-4,7,10,13-tetra-oxa-pentadecanoyl]-L-lysine CAS-NO: 1334172-64-3 FORMULA: C ₄₂ H ₅₉ N ₅ O ₁₁ S MOL WEIGHT: 842,01 g/mol		PEG4440.0100	100 mg	€ 325,00
		PEG4440.1000	1 g	€ 1.300,00
PEG4450 Fmoc-L-Lys(dPEG™(12)-Biotin)-OH N-alpha-(9-Fluorenylmethoxycarbonyl)-N-epsilon-[alpha-Biotin-omega-propionyl dodeca(ethylene glycol)]-L-lysine CAS-NO: 1334172-65-4 FORMULA: C ₅₈ H ₉₁ N ₅ O ₁₉ S MOL WEIGHT: 1194,43 g/mol		PEG4450.0100	100 mg	€ 440,00
		PEG4450.1000	1 g	€ 1.750,00
PEG2166 Biotin-PEG(4)-Lys(PEG(4)-Biotin)-PEG(4)-TFP N-alpha-N-epsilon-Bis(Biotinyl-PEG(4))-lysyl-PEG(4)-2,3,5,6-tetrafluorophenyl ester FORMULA: C ₆₅ H ₁₀₅ F ₄ N ₉ O ₂₁ S ₂ MOL WEIGHT: 1488,7 g/mol		PEG2166.0100	100 mg	€ 275,00
		PEG2166.0001	1 g	€ 1.450,00

Explanation see p.9 (Amine Reactive Biotinylating Reagents)

Desired biotinylating reagent not found in this booklet?

Please ask for our custom synthesis service!



CUSTOM
SYNTHESIS

Fast & efficient from mg to multi-kg
Route Scouting - Process Development - Production


EMPOWERING PEPTIDE INNOVATION

PEG Containing Biotinylation Reagents



PEGs are molecules comprised of repeating ethylene glycol units and are built up in a stepwise fashion. These hydrophilic, non-immunogenic reagents have shown promise in a number of potentially important applications. PEGs have a dramatic impact on protein solubility, stability, and aggregation, and are able to improve ELISAs, IHC assays, and formulations in many aspects. Some protein biotin conjugates are particularly susceptible to aggregation and loss of antigen binding ability, which can complicate assays and provide irreproducible and/or erroneous results.

In comparison, biotinylated IgGs employing a PEG linker show no tendency for aggregation, maintain a high affinity for antigen, and the resulting complexes with the biotin-binding proteins remain completely soluble.

A ...the OLD way!



...the NEW dPEG® way!



Comparison between Biotin-LC labeled and Biotin-dPEG® labeled antibodies:











- > (A) Schematic representation of biotin labeled antibody: "Classical" aliphatic LC linker creates hydrophobic surface with biotin moiety in very close proximity (left side of IgG scheme). dPEG®12 linker creates hydrophilic "slippery" surface with biotin moiety given more space for interaction (right side of IgG scheme).
- > (B) Mixing 3 nmol streptavidin with 0.9 nmol either Biotin-LC-GAR antibody (left tube) or Biotin-dPEG®12-GAR antibody (right tube) clearly demonstrates the ability of the hydrophilic dPEG® spacer to solubilize the resulting streptavidin-biotin-IgG complex completely while the streptavidin-biotin-LC-GAR complex is obviously precipitating from solution.



Biotin-PEG-Acids / Amine Reactive Biotinylation Reagents

Active esters such as N-hydroxysuccinimide (NHS) esters permit a catalyst-free conjugation reaction to amines under mild conditions. Yet, NHS esters may hydrolyze over time if, for example, atmospheric moisture has condensed on the substance.

As alternative, tetrafluorophenyl (TFP) and pentafluorophenyl (PFP) active esters can be applied that are much more resistant to spontaneous hydrolysis.

	Article No.	Quantity	Price
PEG2560 Biotin-O₂Oc-OH 8-Biotinylamido-3,6-dioxaoctanoic acid CAS-NO: 1238575-77-3 FORMULA: C ₁₆ H ₂₇ N ₃ O ₆ S MOL WEIGHT: 389,47 g/mol	PEG2560.0250	250 mg	€ 100,00
	PEG2560.0500	500 mg	€ 180,00
	PEG2560.0001	1 g	€ 280,00
	PEG2560.0005	5 g	€ 1.000,00
PEG2565 Biotin-O₂Oc-O₂Oc-OH 8-Biotinylamido-3,6-dioxaoctanoic acid dimer FORMULA: C ₂₂ H ₃₈ N ₄ O ₈ S MOL WEIGHT: 534,62 g/mol	PEG2565.0100	100 mg	€ 120,00
	PEG2565.0250	250 mg	€ 200,00
	PEG2565.0500	500 mg	€ 360,00
	PEG2565.0001	1 g	€ 560,00
	PEG2565.0005	5 g	€ 2.000,00
PEG5280 Biotin-TOTA-glutaramic acid*DIPEA 5,21-dioxo-25-((3aS,4S,6aR)-2-oxohexahydro-1H-thieno[3,4-d]imidazol-4-yl)-10,13,16-trioxo-6,20-diazapentacosan-1-oic acid diisopropylethylamine salt FORMULA: C ₂₅ H ₄₄ N ₄ O ₈ S*C ₈ H ₁₇ N MOL WEIGHT: 560,70*129,30 g/mol	PEG5280.0100	100 mg	€ 96,00
	PEG5280.0250	250 mg	€ 160,00
	PEG5280.0500	500 mg	€ 288,00
	PEG5280.1000	1 g	€ 448,00
	PEG5280.5000	5 g	€ 1.600,00

		Article No.	Quantity	Price
PEG2071	Biotin-TEG-TFP	PEG2071.0100	100 mg	€ 225,00
		PEG2071.1000	1 g	€ 825,00
Biotin-tetra(ethylene glycol)-2,3,5,6-tetrafluorophenyl ester FORMULA: $C_{27}H_{37}F_4N_3O_8S$ MOL WEIGHT: 639,66 g/mol				
PEG1515	Biotin-PEG(4)-COOH	PEG1515.0250	250 mg	€ 225,00
		PEG1515.0001	1 g	€ 675,00
15-Biotinamino-4,7,10,13-tetraoxa-pentadecanoic acid CAS-NO: 721431-18-1 FORMULA: $C_{21}H_{33}N_3O_5S$ MOL WEIGHT: 491,6 g/mol				
PEG1870	Biotin-dPEG(4)-NHS	PEG1870.0050	50 mg	€ 225,00
		PEG1870.0001	1 g	€ 990,00
15-Biotinamino-4,7,10,13-tetraoxa-pentadecanoic acid succinimidyl ester CAS-NO: 459426-22-3 FORMULA: $C_{25}H_{40}N_4O_{10}S$ MOL WEIGHT: 588,67 g/mol				
PEG1845	Biotin-dPEG(4)-NHS-(Biotinidase resistant)	PEG1845.0100	100 mg	€ 275,00
		PEG1845.0001	1 g	€ 1.550,00
18-Biotinamino-17-oxo-4,7,10,13-tetraoxa-16-azaiicosan-1-oic acid succinimidyl ester CAS-NO: 1334172-61-0 FORMULA: $C_{29}H_{47}N_5O_{11}S$ MOL WEIGHT: 673,78 g/mol				
PEG7740	Biotin-dPEG™(4)-TFP	PEG7740.0100	100 mg	€ 200,00
		PEG7740.0001	1 g	€ 825,00
alpha-Biotin-omega-(2,3,5,6-tetrafluorophenyl propionate) tetra(ethylene glycol) FORMULA: $C_{27}H_{37}F_4N_3O_8S$ MOL WEIGHT: 639,66 g/mol				
PEG7760	Biotin-dPEG™(4)-TFP-(Biotinidase resistant)	PEG7760.0100	100 mg	€ 275,00
		PEG7760.0001	1 g	€ 1.550,00
18-Biotinamino-17-oxo-4,7,10,13-tetraoxa-16-azaiicosan-1-oic acid (2,3,5,6-tetrafluorophenyl) ester FORMULA: $C_{31}H_{47}F_4N_5O_{11}S$ MOL WEIGHT: 724,76 g/mol				
PEG6690	Biotin-PEG(6)-COOH	PEG6690.0000	please inquire	
21-[D(+)-Biotinylamino]-4,7,10,13,16,19-hexaoxaheneicosanoic acid CAS-NO: 1352814-10-8 FORMULA: $C_{25}H_{45}N_3O_{10}S$ MOL WEIGHT: 579,7 g/mol				
PEG1051	Biotin-PEG(12)-COOH	PEG1051.0100	100 mg	€ 290,00
		PEG1051.0001	1 g	€ 1.250,00
alpha-Biotin-omega-(propionic acid)-dodeca(ethylene glycol) CAS-NO: 948595-11-7 FORMULA: $C_{37}H_{69}N_3O_{16}S$ MOL WEIGHT: 844,04 g/mol				
PEG1860	Biotin-dPEG(12)-NHS	PEG1860.0050	50 mg	€ 275,00
		PEG1860.0001	1 g	€ 1.600,00
alpha-Biotin-omega-carboxy succinimidyl ester dodeca(ethylene glycol) CAS-NO: 365441-71-0 FORMULA: $C_{41}H_{72}N_4O_{18}S$ MOL WEIGHT: 941,09 g/mol				
PEG5130	Biotin-PEG(12)-TFP	PEG5130.0100	100 mg	€ 275,00
		PEG5130.1000	1 g	€ 1.275,00
alpha-Biotin-omega-propionic acid (2,3,5,6-tetrafluorophenyl) ester dodeca(ethylene glycol) FORMULA: $C_{43}H_{69}F_4N_3O_{16}S$ MOL WEIGHT: 992,08 g/mol				












		Article No.	Quantity	Price
PEG7770	Biotin-dPEG™(12)-TFP-(Biotinidase resistant)	PEG7770.0100	100 mg	€ 380,00
alpha-Biotinyl-(2-aminobutyl)-omega-(2,3,5,6-tetrafluorophenyl propionate) dodeca (ethylene glycol) FORMULA: $C_{47}H_{76}F_4N_4O_{17}S$ MOL WEIGHT: 1077,18 g/mol		PEG7770.0001	1 g	€ 2.300,00
				
PEG4260	Biotin-dPEG™(24)-COOH	PEG4260.0100	100 mg	€ 350,00
alpha-Biotin-omega-(propionic acid) 24(ethylene glycol) CAS-NO: 721431-18-1 FORMULA: $C_{61}H_{117}N_3O_{28}S$ MOL WEIGHT: 1372.65 g/mol		PEG4260.1000	1 g	€ 1.750,00
				
PEG4250	Biotin-dPEG™(24)-NHS	PEG4250.0050	50 mg	€ 380,00
alpha-Biotin-omega-(succinimidyl propionate) 24(ethylene glycol) CAS-NO: 365441-71-0 FORMULA: $C_{65}H_{120}N_4O_{30}S$ MOL WEIGHT: 1469,72 g/mol		PEG4250.1000	1 g	€ 1.800,00
				
PEG7750	Biotin-dPEG™(24)-TFP	PEG7750.0100	100 mg	€ 380,00
alpha-Biotin-omega-(2,3,5,6-tetrafluorophenyl propionate) 24(ethylene glycol) FORMULA: $C_{67}H_{117}F_4N_3O_{28}S$ MOL WEIGHT: 1520,71 g/mol		PEG7750.0001	1 g	€ 1.800,00
				
PEG4270	Biotin-dPEG™(48)-COOH	PEG4270.0100	100 mg	€ 440,00
alpha-Biotin-omega-(propionic acid) 48(ethylene glycol) CAS-NO: 721431-18-1 FORMULA: $C_{112}H_{218}N_4O_{53}S$ MOL WEIGHT: 2500,99 g/mol		PEG4270.1000	1 g	€ 1.950,00
				
PEG1053	Biotin-PEG-COOH	PEG1053.0500	500 mg	€ 450,00
alpha-Biotin-omega-carboxy poly(ethylene glycol) (PEG-MW 3.000 Da) MOL WEIGHT: 3000 Da		PEG1053.0001	1 g	€ 790,00
PEG1054	Biotin-PEG-COOH	PEG1054.0500	500 mg	€ 450,00
alpha-Biotin-omega-carboxy poly(ethylene glycol) (PEG-MW 5.000 Da) MOL WEIGHT: 5000 Da		PEG1054.0001	1 g	€ 790,00
				
PEG1052	Biotin-PEG-COOH	PEG1052.0500	500 mg	€ 480,00
alpha-Biotin-omega-carboxy poly(ethylene glycol) (PEG-MW 10.000 Da) MOL WEIGHT: 10000 Da		PEG1052.0001	1 g	€ 850,00
PEG1056	Biotin-PEG-NHS	PEG1056.0100	100 mg	€ 140,00
alpha-Biotin-omega-carboxy succinimidyl ester poly(ethylene glycol) (PEG-MW 3.000 Da) MOL WEIGHT: 3000 Da		PEG1056.0500	500 mg	€ 510,00
		PEG1056.0001	1 g	€ 910,00
PEG1057	Biotin-PEG-NHS	PEG1057.0500	500 mg	€ 510,00
alpha-Biotin-omega-carboxy succinimidyl ester poly(ethylene glycol) (PEG-MW 5.000 Da) MOL WEIGHT: 5000 Da		PEG1057.0001	1 g	€ 910,00
				
PEG1055	Biotin-PEG-NHS	PEG1055.0500	500 mg	€ 550,00
alpha-Biotin-omega-carboxy succinimidyl ester poly(ethylene glycol) (PEG-MW 10.000 Da) MOL WEIGHT: 10000 Da		PEG1055.0001	1 g	€ 970,00

		Article No.	Quantity	Price
PEG1580 mal-dPEG(4)-NHHH-Boc	1-Maleimido-3-oxo-7,10,13,16-tetraoxa-4-azanadecan-19-oyl-(N'-t-butylloxycarbonyl) hydrazid CAS-NO: 1127247-28-2 FORMULA: C ₂₃ H ₃₈ N ₄ O ₁₀ MOL WEIGHT: 530,57 g/mol			
PEG1590 mal-dPEG(8)-NHS	1-Maleimido-3-oxo-7,10,13,16,19,22,25,28-octaoxa-4-azahentriacontan-31-oic acid succinimidyl ester CAS-NO: 756525-93-6 FORMULA: C ₃₀ H ₄₇ N ₃ O ₁₅ MOL WEIGHT: 689,71 g/mol			
PEG1600 mal-dPEG(24)-COOH	alpha-Maleimido-24(ethylene glycol)-omega-propionic acid CAS-NO: 871133-36-7 FORMULA: C ₅₈ H ₁₀₈ N ₂ O ₂₉ MOL WEIGHT: 1297,47 g/mol			
PEG6590 Biotin-PEG-SQA	alpha-Biotin-omega-squaric acid ethyl ester poly(ethylene glycol) (PEG-MW 3000 Dalton) MOL WEIGHT: 3000 Da			
PEG6595 Biotin-PEG-SQA	alpha-Biotin-omega-squaric acid ethyl ester poly(ethylene glycol) (PEG-MW 5000 Dalton) MOL WEIGHT: 5000 Da			
PEG6580 Biotin-PEG-SQA	alpha-Biotin-omega-squaric acid ethyl ester poly(ethylene glycol) (PEG-MW 10000 Dalton) MOL WEIGHT: 10000 Da			
PEG6585 Biotin-PEG-SQA	alpha-Biotin-omega-squaric acid ethyl ester poly(ethylene glycol) (PEG-MW 20000 Dalton) MOL WEIGHT: 20000 Da			





Squaric acid ethyl esters for chemoselective biotinylation of lysine side chains in proteins.

Biotin-PEG-Amines

		Article No.	Quantity	Price
PEG2760 Biotin-DOOA*HCl	1-Biotinyl-3,6-dioxa-8-octaneamine hydrochloride CAS-NO: 862373-14-6 FORMULA: C ₁₆ H ₃₀ N ₄ O ₅ S*HCl MOL WEIGHT: 374,50*36,45 g/mol			
PEG2110 Biotin-TEG-NH₂*TFA	N-Biotin-tetra(ethylene glycol)-diamine CAS-NO: 1334172-59-6 FORMULA: C ₂₀ H ₃₈ N ₄ O ₅ S*CF ₃ CO ₂ H MOL WEIGHT: 446,60*114,02 g/mol			
PEG5140 Biotin-TEG-O-NH₂*HCl	N-(2-(2-(2-(2-(aminoxy)ethoxy)ethoxy)ethoxy)ethyl)-5-biotinamide hydrochloride FORMULA: C ₁₈ H ₃₄ N ₄ O ₆ S*HCl MOL WEIGHT: 434,55*36,45 g/mol			

		Article No.	Quantity	Price
PEG7780	Biotin-dPEG™(3)-NH₂*TFA	PEG7780.0100	100 mg	€ 200,00
		PEG7780.0001	1 g	€ 990,00
<p>alpha-Biotin-omega-amino-tri(ethylene glycol) trifluoroacetate</p> <p>FORMULA: C₂₀H₃₈N₄O₅S*CF₃CO₂H</p> <p>MOL WEIGHT: 446,60*114,01 g/mol</p>				
PEG7790	Biotin-dPEG™(3)-O-NH₂*HCl	PEG7790.0050	50 mg	€ 200,00
		PEG7790.0001	1 g	€ 990,00
<p>FORMULA: C₁₈H₃₄N₄O₆*HCl</p> <p>MOL WEIGHT: 434,55*36,45 g/mol</p>				
PEG1425	Biotin-dPEG(4)-NHNH₂	PEG1425.0050	50 mg	€ 380,00
		PEG1425.0001	1 g	€ 1.450,00
<p>15-Biotinamino-4,7,10,13-tetraoxa-pentadecanoyl hydrazide</p> <p>CAS-NO: 756525-97-0</p> <p>FORMULA: C₂₁H₃₉N₅O₅</p> <p>MOL WEIGHT: 505,63 g/mol</p>				
PEG4290	Biotin-dPEG™(7)-NH₂	PEG4290.0100	100 mg	€ 380,00
		PEG4290.1000	1 g	€ 1.200,00
<p>alpha-Biotin-omega-amino hepta(ethylene glycol)</p> <p>CAS-NO: 1334172-76-7</p> <p>FORMULA: C₂₆H₅₀N₄O₅S</p> <p>MOL WEIGHT: 594,76 g/mol</p>				
PEG1044	Biotin-PEG(9)-NH₂	PEG1044.0100	100 mg	€ 390,00
		PEG1044.0001	1 g	€ 1.160,00
<p>alpha-Biotin-omega-amino-nona(ethylene glycol)</p> <p>FORMULA: C₃₀H₅₈N₄O₁₁S</p> <p>MOL WEIGHT: 682,88 g/mol</p>				
PEG4300	Biotin-dPEG™(11)-NH₂	PEG4300.0100	100 mg	€ 425,00
		PEG4300.1000	1 g	€ 1.300,00
<p>alpha-Biotin-omega-amino undeca(ethylene glycol)</p> <p>CAS-NO: 604786-74-5</p> <p>FORMULA: C₃₄H₆₆N₄O₁₃S</p> <p>MOL WEIGHT: 770,97 g/mol</p>				
PEG4680	Biotin-dPEG™(11)-O-NH₂*HCl	PEG4680.0050	50 mg	€ 275,00
		PEG4680.1000	1 g	€ 1.600,00
<p>alpha-Biotinyl-omega-oxyamine-undeca(ethylene glycol) hydrochloride</p> <p>FORMULA: C₃₄H₆₆N₄O₁₄S</p> <p>MOL WEIGHT: 786,97*36,45 g/mol</p>				
PEG4310	Biotin-dPEG™(23)-NH₂	PEG4310.0100	100 mg	€ 440,00
		PEG4310.1000	1 g	€ 1.450,00
<p>alpha-Biotin-omega-amino 23(ethylene glycol)</p> <p>CAS-NO: 604786-74-5</p> <p>FORMULA: C₅₈H₁₁₄N₄O₂₅S</p> <p>MOL WEIGHT: 1299,6 g/mol</p>				
PEG1046	Biotin-PEG-NH₂	PEG1046.0500	500 mg	€ 410,00
		PEG1046.0001	1 g	€ 670,00
<p>alpha-Biotin-omega-amino poly(ethylene glycol) (PEG-MW 3.000 Dalton)</p> <p>MOL WEIGHT: 3000 Da</p>				
PEG1047	Biotin-PEG-NH₂	PEG1047.0500	500 mg	€ 410,00
		PEG1047.0001	1 g	€ 670,00
<p>alpha-Biotin-omega-amino poly(ethylene glycol) (PEG-MW 5.000 Dalton)</p> <p>MOL WEIGHT: 5000 Da</p>				
PEG1045	Biotin-PEG-NH₂	PEG1045.0500	500 mg	€ 440,00
		PEG1045.0001	1 g	€ 730,00
<p>alpha-Biotin-omega-amino poly(ethylene glycol) (PEG-MW 10.000 Dalton)</p> <p>MOL WEIGHT: 10000 Da</p>				

Biotin-PEG-Maleimides

		Article No.	Quantity	Price
PEG1605 Biotin-dPEG(3)-mal N-(3-(2-(2-(3-(Biotinamino)propoxy)ethoxy)ethoxy)propyl)-3-maleimidylpropanamide CAS-NO: 525573-22-2 FORMULA: C ₂₇ H ₄₃ N ₅ O ₈ S MOL WEIGHT: 597,72 g/mol		PEG1605.0250	250 mg	€ 225,00
		PEG1605.0001	1 g	€ 650,00
PEG1595 Biotin-dPEG(11)-mal alpha-Biotin-omega-maleimido undeca(ethylene glycol) CAS-NO: 1334172-60-9 FORMULA: C ₄₁ H ₇₇ N ₅ O ₁₀ S MOL WEIGHT: 922,09 g/mol		PEG1595.0025	25 mg	€ 350,00
		PEG1595.0100	100 mg	€ 625,00
PEG4320 Biotin-dPEG™(23)-mal alpha-Biotin-omega-maleimido 23(ethylene glycol) CAS-NO: 1334172-60-9 FORMULA: C ₆₅ H ₁₁₉ N ₅ O ₂₈ S MOL WEIGHT: 1450,72 g/mol		PEG4320.0025	25 mg	€ 380,00
		PEG4320.0100	100 mg	€ 690,00
PEG1049 Biotin-PEG-mal alpha-Biotin-omega-maleimido poly(ethylene glycol) (PEG-MW 3.000 Dalton) MOL WEIGHT: 3000 Da		PEG1049.0500	500 mg	€ 550,00
		PEG1049.0001	1 g	€ 910,00
PEG1050 Biotin-PEG-mal alpha-Biotin-omega-maleimido poly(ethylene glycol) (PEG-MW 5.000 Dalton) MOL WEIGHT: 5000 Da		PEG1050.0500	500 mg	€ 550,00
		PEG1050.0001	1 g	€ 910,00
PEG1048 Biotin-PEG-mal alpha-Biotin-omega-maleimido poly(ethylene glycol) (PEG-MW 10.000 Dalton) MOL WEIGHT: 10000 Da		PEG1048.0500	500 mg	€ 610,00
		PEG1048.0001	1 g	€ 970,00

Interested in further PEGylating reagents?







Find many more in our brochure "Comprehensive Drug Delivery Survey"!

Iris
BIOTECH GMBH

COMPREHENSIVE
DRUG DELIVERY SURVEY









FROM GRAMS TO MULTI-TON LOTS

Biotin-PEG-Thiols

		Article No.	Quantity	Price
PEG3570 Lipoamide-dPEG™(3)-biotin alpha-Lipoamide-omega-biotinyl tri(ethylene glycol) CAS-NO: 1334172-74-5 FORMULA: C ₂₈ H ₅₀ N ₄ O ₆ S ₃ MOL WEIGHT: 634,91 g/mol		PEG3570.0100	100 mg	€ 380,00
		PEG3570.1000	1 g	€ 1.275,00
PEG4705 Biotin-PEG(11)-SH alpha-Biotin-omega-mercapto undeca(ethylene glycol) FORMULA: C ₃₄ H ₆₅ N ₃ O ₁₃ S ₂ MOL WEIGHT: 788,02 g/mol		PEG4705.0100	100 mg	€ 450,00
		PEG4705.1000	1 g	€ 1.380,00
PEG3580 Lipoamide-dPEG™(11)-biotin alpha-Lipoamide-omega-biotinyl undeca(ethylene glycol) CAS-NO: 960069-81-2 FORMULA: C ₄₂ H ₇₈ N ₄ O ₁₄ S ₃ MOL WEIGHT: 959,28 g/mol		PEG3580.0100	100 mg	€ 500,00
		PEG3580.1000	1 g	€ 1.600,00
PEG7580 mal-dPEG™(24)-TFP alpha-Maleinimido-24(ethylen glycol)-omega-propionic acid 2,3,5,6-tetrafluorophenyl ester FORMULA: C ₆₄ H ₁₀₈ F ₄ N ₂ O ₂₉ MOL WEIGHT: 1445,53 g/mol		PEG7580.0100	100 mg	€ 380,00
		PEG7580.0001	1 g	€ 1.600,00
PEG1213 Biotin-PEG-SH alpha-Biotin-omega-mercapto poly(ethylene glycol) (PEG-MW 3.000 Dalton) MOL WEIGHT: 3000 Da		PEG1213.0500	500 mg	€ 900,00
		PEG1213.0001	1 g	€ 1.580,00
PEG1214 Biotin-PEG-SH alpha-Biotin-omega-mercapto poly(ethylene glycol) (PEG-MW 5.000 Dalton) MOL WEIGHT: 5000 Da		PEG1214.0500	500 mg	€ 900,00
		PEG1214.0001	1 g	€ 1.580,00
PEG1212 Biotin-PEG-SH alpha-Biotin-omega-mercapto poly(ethylene glycol) (PEG-MW 10.000 Dalton) MOL WEIGHT: 10000 Da		PEG1212.0500	500 mg	€ 930,00
		PEG1212.0001	1 g	€ 1.640,00
PEG1226 Biotin-PEG-SH alpha-Biotin-omega-mercapto poly(ethylene glycol) (PEG-MW 20.000 Dalton) MOL WEIGHT: 20000 Da		PEG1226.0000	please inquire	
PEG4770 Biotin-PEG-OPSS alpha-Biotin-omega-pyridyl-2-disulfid poly(ethylene glycol) (PEG-MW 3.000 Dalton) MOL WEIGHT: 3000 Da		PEG4770.0500	500 mg	€ 900,00
		PEG4770.1000	1 g	€ 1.580,00
PEG4775 Biotin-PEG-OPSS alpha-Biotin-omega-pyridyl-2-disulfid poly(ethylene glycol) (PEG-MW 5.000 Dalton) MOL WEIGHT: 5000 Da		PEG4775.0500	500 mg	€ 900,00
		PEG4775.1000	1 g	€ 1.580,00
PEG4780 Biotin-PEG-OPSS alpha-Biotin-omega-pyridyl-2-disulfid poly(ethylene glycol) (PEG-MW 10.000 Dalton) MOL WEIGHT: 10000 Da		PEG4780.0500	500 mg	€ 900,00
		PEG4780.1000	1 g	€ 1.580,00
PEG4785 Biotin-PEG-OPSS alpha-Biotin-omega-pyridyl-2-disulfid poly(ethylene glycol) (PEG-MW 20.000 Dalton) MOL WEIGHT: 20000 Da		PEG4785.0500	500 mg	€ 900,00
		PEG4785.1000	1 g	€ 1.580,00

PEG OPSS reagents for reversible biotinylation of thiols.

Other Biotin-PEG Reagents

		Article No.	Quantity	Price
PEG1420	Biotin-dPEG™(3)-Cyanocobalamin	PEG1420.0005	5 mg	€ 250,00
<p>N-(3-(2-(2-(3-(Biotinamino)propoxy)ethoxy)ethoxy)propyl)-cyanocobalamin</p> <p>CAS-NO: 295329-79-2 FORMULA: C₈₃H₁₂₃N₁₇O₂₀CoPS MOL WEIGHT: 1800,94 g/mol</p>				
PEG2080	Biotin-TEG-Biotin	PEG2080.0050	50 mg	€ 225,00
<p>N,N'-Bisbiotin-tetra(ethylene glycol)-diamine</p> <p>CAS-NO: 194920-54-2 FORMULA: C₃₀H₅₂N₆O₈S₂ MOL WEIGHT: 672,9 g/mol</p>				
PEG1206	Biotin-PEG-OH	PEG1206.0500	500 mg	€ 720,00
<p>alpha-Biotin-omega-hydroxy poly(ethylene glycol) (PEG-MW 3.000 Dalton)</p> <p>MOL WEIGHT: 3000 Da</p>				
PEG1207	Biotin-PEG-OH	PEG1207.0500	500 mg	€ 720,00
<p>alpha-Biotin-omega-hydroxy poly(ethylene glycol) (PEG-MW 5.000 Dalton)</p> <p>MOL WEIGHT: 5000 Da</p>				
PEG4850	Biotin-PEG-Si(OMe)3	PEG4850.0500	500 mg	€ 560,00
<p>alpha-Biotinyl-omega-trimethoxysilyl poly(ethylene glycol) (PEG-MW 3.000 Dalton)</p> <p>MOL WEIGHT: 3000 Da</p>				
PEG4855	Biotin-PEG-Si(OMe)3	PEG4855.0500	500 mg	€ 560,00
<p>alpha-Biotinyl-omega-trimethoxysilyl poly(ethylene glycol) (PEG-MW 5.000 Dalton)</p> <p>MOL WEIGHT: 5000 Da</p>				
PEG4860	Biotin-PEG-Si(OMe)3	PEG4860.0500	500 mg	€ 560,00
<p>alpha-Biotinyl-omega-trimethoxysilyl poly(ethylene glycol) (PEG-MW 10.000 Dalton)</p> <p>MOL WEIGHT: 10000 Da</p>				
PEG4865	Biotin-PEG-Si(OMe)3	PEG4865.0500	500 mg	€ 560,00
<p>alpha-Biotinyl-omega-trimethoxysilyl poly(ethylene glycol) (PEG-MW 20.000 Dalton)</p> <p>MOL WEIGHT: 20000 Da</p>				
<p>PEG silanes for biotinylation of silicate surfaces.</p>				

2. Terms and Conditions of Sales

All orders placed by a buyer are accepted and all contracts are made subject to the terms which shall prevail and be effective notwithstanding any variations or additions contained in any order or other document submitted by the buyer. No modification of these terms shall be binding upon Iris Biotech GmbH unless made in writing by an authorised representative of Iris Biotech GmbH.

Placing of Orders

Every order made by the buyer shall be deemed an offer by the buyer to purchase products from Iris Biotech GmbH and will not be binding on Iris Biotech GmbH until a duly authorised representative of Iris Biotech GmbH has accepted the offer made by the buyer. Iris Biotech GmbH may accept orders from commercial, educational or government organisations, but not from private individuals and Iris Biotech GmbH reserves the right to insist on a written order and/or references from the buyer before proceeding.

There is no minimum order value. At the time of acceptance of an order Iris Biotech GmbH will either arrange prompt despatch from stock or the manufacture/acquisition of material to satisfy the order. In the event of the latter Iris Biotech GmbH will indicate an estimated delivery date. In addition to all its other rights Iris Biotech GmbH reserves the right to refuse the subsequent cancellation of the order if Iris Biotech GmbH expects to deliver the product on or prior to the estimated delivery date. Time shall not be of the essence in respect of delivery of the products. If Iris Biotech GmbH is unable to deliver any products by reason of any circumstances beyond its reasonable control („Force Majeure“) then the period for delivery shall be extended by the time lost due to such Force Majeure. Details of Force Majeure will be forwarded by Iris Biotech GmbH to the buyer as soon as reasonably practicable.

Prices, Quotations and Payments

Prices are subject to change. For the avoidance of doubt, the price advised by Iris Biotech GmbH at the time of the buyer placing the order shall supersede any previous price indications. The buyer must contact the local office of Iris Biotech GmbH before ordering if further information is required. Unless otherwise agreed by the buyer and Iris Biotech GmbH, the price shall be for delivery ex-works. In the event that the buyer requires delivery of the products otherwise than ex-works the buyer should contact the local office of Iris Biotech GmbH in order to detail its requirements. Iris Biotech GmbH shall, at its discretion, arrange the buyer's delivery requirements including, without limitation, transit insurance, the mode of transit (Iris Biotech GmbH reserves the right to vary the mode of transit if any regulations or other relevant considerations so require) and any special packaging requirements (including cylinders). For the avoidance of doubt all costs of delivery and packaging in accordance with the buyer's requests over and above that of delivery in standard packaging ex-works shall be for the buyer's account unless otherwise agreed by both parties. Incoterms 2010 shall apply. Any tax, duty or charge imposed by governmental authority or otherwise and any other applicable taxes, duties or charges shall be for the buyer's account. Iris Biotech GmbH may, on request and where possible, provide quotations for multiple packs or bulk quantities, and non-listed items. Irrespective of the type of request or means of response all quotations must be accepted by the buyer without condition and in writing before an order will be accepted by Iris Biotech GmbH. Unless agreed in writing on different terms, quotations are valid for 30 days from the date thereof. Payment terms are net 30 days from invoice date unless otherwise agreed in writing. Iris Biotech GmbH reserves the right to request advance payment at its discretion. For overseas transactions the buyer shall pay all the banking charges of Iris Biotech GmbH. The buyer shall not be entitled to withhold or set-off payment for the products for any reason whatsoever. Government/Corporate Visa and MasterCard (and other such credit cards) may be accepted on approved accounts for payment of the products. Personal credit cards are not acceptable. Failure to comply with the terms of payment of Iris Biotech GmbH shall constitute default without reminder. In these circumstances Iris Biotech GmbH may (without prejudice to any other of its rights under these terms) charge interest to accrue on a daily basis at the rate of 2% per month from the date upon which payment falls due to the actual date of payment (such interest shall be paid monthly). If the buyer shall fail to fulfil the payment terms in respect of any invoice of Iris Biotech GmbH Iris Biotech GmbH may demand payment of all outstanding balances from the buyer whether due or not and/or cancel all outstanding orders and/or decline to make further deliveries or provision of services except upon receipt of cash or satisfactory securities. Until payment by the buyer in full of the price and any other monies due to Iris Biotech GmbH in respect of all other products or services supplied or agreed to be supplied by Iris Biotech GmbH to the buyer (including but without limitation any costs of delivery) the property in the products shall remain vested in Iris Biotech GmbH.

Shipping, Packaging and Returns

The buyer shall inspect goods immediately on receipt and inform Iris Biotech GmbH of any shortage or damage within five days. Quality problems must be notified within ten days of receipt. Goods must not be returned without prior written authorisation of Iris Biotech GmbH. Iris Biotech GmbH shall at its sole discretion replace the defective products (or parts thereof) free of charge or refund the price (or proportionate price) to buyer. Opened or damaged containers cannot be returned by the buyer without the written prior agreement of Iris Biotech GmbH. In the case of agreed damaged containers which cannot be so returned, the buyer assumes responsibility for the safe disposal of such containers in accordance with all applicable laws.

Product Quality, Specifications and Technical Information

Products are analysed in the Quality Control laboratories of Iris Biotech GmbH's production partners by methods and procedures which Iris Biotech GmbH considers appropriate. In the event of any dispute concerning reported discrepancies arising from the buyer's analytical results, determined by the buyer's own analytical procedures, Iris Biotech GmbH reserves the right to rely on the results of own analytical methods of Iris Biotech GmbH. Certificates of Analysis or Certificates of Conformity are available at the discretion of Iris Biotech GmbH for bulk orders but not normally for prepack orders. Iris Biotech GmbH reserves the right to make a charge for such Certification. Specifications may change and reasonable variation from any value listed should not form the basis of a dispute. Any supply by Iris Biotech GmbH of bespoke or custom product for a buyer shall be to a specification agreed by both parties in writing. Technical information, provided orally, in writing, or by electronic means by or on behalf of Iris Biotech GmbH, including any descriptions, references, illustrations or diagrams in any Catalogue or brochure, is provided for guidance purposes only and is subject to change.

Safety

All chemicals should be handled only by competent, suitably trained persons, familiar with laboratory procedures and potential chemical hazards. The burden of safe use of the products of Iris Biotech GmbH vests in the buyer. The buyer assumes all responsibility for warning his employees, and any persons who might reasonably be expected to come into contact with the products, of all risks to person and property in any way connected with the products and for instructing them in their safe handling and use. The buyer also assumes the responsibility for the safe disposal of all products in accordance with all applicable laws.

Uses, Warranties and Liabilities

All products of Iris Biotech GmbH are intended for laboratory research purposes and unless otherwise stated on product labels, in the catalogue and product information sheet of Iris Biotech GmbH or in other literature furnished to the buyer, are not to be used for any other purposes, including but not limited to use as or as components in drugs for human or animal use, medical devices, cosmetics, food additives, household chemicals, agricultural or horticultural products or pesticides. Iris Biotech GmbH offers no warranty regarding the fitness of any product for a particular purpose and shall not be responsible for any loss or damage whatsoever arising there from. No warranty or representation is given by Iris Biotech GmbH that the products do not infringe any letters patent, trademarks, registered designs or other industrial rights. The buyer further warrants to Iris Biotech GmbH that any use of the products in the United States of America shall not result in the products becoming adulterated or misbranded within the meaning of the Federal Food, Drug and Cosmetic Act (or such equivalent legislation in force in the buyer's jurisdiction) and shall not be materials which may not, under sections 404, 505 or 512 of the Act, be introduced into interstate commerce. The buyer acknowledges that, since the products of Iris Biotech GmbH are intended for research purposes, they may not be on the Toxic Substances Control Act 1976 („TSCA“) inventory. The buyer warrants that it shall ensure that the products are approved for use under the TSCA (or such other equivalent legislation in force in the buyer's jurisdiction), if applicable. The buyer shall be responsible for complying with any legislation or regulations governing the use of the products and their importation into the country of destination (for the avoidance of doubt to include, without limitation, the TSCA and all its amendments, all EINECS, ELINCS and NONS regulations). If any licence or consent of any government or other authority shall be required for the acquisition, carriage or use of the products by the buyer the buyer shall obtain the same at its own expense and if necessary produce evidence of the same to Iris Biotech GmbH on demand. Failure to do so shall not entitle the buyer to withhold or delay payment. Any additional expenses or charges incurred by Iris Biotech GmbH resulting from such failure shall be for the buyer's account. Save for death or personal injury caused by negligence of Iris Biotech GmbH, sole obligation of Iris Biotech GmbH and buyer's exclusive remedy with respect to the products proved to the satisfaction of Iris Biotech GmbH to be defective or products incorrectly supplied shall be to accept the return of said products to Iris Biotech GmbH for refund of the actual purchase price paid by the buyer (or proportionate part thereof), or replacement of the defective product (or part thereof) with alternative product. Iris Biotech GmbH shall have no liability to the buyer under or arising directly or indirectly out of or otherwise in connection with the supply of products by Iris Biotech GmbH to the buyer and/or their re-sale or use by the buyer or for any product, process or services of the buyer which in any way comprises the product in contract tort (including negligence or breach of statutory duty) or otherwise for pure economic loss, loss of profit, business, reputation, depletion of brand, contracts, revenues or anticipated savings or for any special indirect or consequential damage or loss of any nature except as may otherwise be expressly provided for in these terms. All implied warranties, terms and representations in respect of the products (whether implied by statute or otherwise) are excluded to the fullest extent permitted by law. The buyer shall indemnify Iris Biotech GmbH for and against any and all losses, damages and expenses, including legal fees and other costs of defending any action, that Iris Biotech GmbH may sustain or incur as a result of any act or omission by the buyer, its officers, agents or employees, its successors or assignees, its customers or all other third parties, whether direct or indirect, in connection with the use of any product. For the avoidance of doubt and in the event that Iris Biotech GmbH supplies bespoke or custom product to the buyer's design or specification, this indemnity shall extend to include any claim by a third party that the manufacture of the product for the buyer or the use of the product by the buyer infringes the intellectual property rights of any third party.

General

Iris Biotech GmbH shall be entitled to assign or sub-contract all or any of its rights and obligations hereunder. The buyer shall not be entitled to assign, transfer, sub-contract or otherwise delegate any of its rights or obligations hereunder. Any delay or forbearance by Iris Biotech GmbH in exercising any right or remedy under these terms shall not constitute a waiver of such right or remedy. If any provision of these terms is held by any competent authority to be invalid or unenforceable in whole or in part the validity of the other provisions of these terms and the remainder of the provision in question shall not be affected. These terms shall be governed by German Law and the German Courts shall have exclusive jurisdiction for the hearing of any dispute between the parties save in relation to enforcement where the jurisdiction of the German Courts shall be non-exclusive.

Publications from Iris Biotech:



We are carrying out many **Contract Manufacturing** projects in the areas of our activities; our strong points are unusual derivatives with one or several chiral centres. Specific propriety platforms give access to a wide variety of unusual amino acid building blocks, polymanies and resins.



Our comprehensive **Guideline on Resins** for solid-supported peptide synthesis. Complete with protocols for resin loading and linker cleavage, standard protocols for peptide chemistry, as well as a current product list of our solid supports. It also includes latest resin developments, e.g. SEA and Hydrazon resins, very useful tools for the synthesis of long peptides via native chemical ligation (NCL).



Brochure for **Click Chemistry** containing over 300 different Azido and Alkyne Compounds for drug discovery, drug delivery and diagnostics.



Disulfide bridges are an important structural element in many proteins and peptides. This brochure is a guideline how to plan and execute the **Synthesis of Cyclic Peptides** with one or several disulfide bridges.



Our **Comprehensive Drug Delivery Survey** contains the latest state-of-the-art collection of carriers for Polymer Therapeutics available today. More than 900 products from the areas PEGylation, poly (amino acids) and PEG-based multifunctional Dendrimers give the medicinal chemist all options to find the right delivery technology from small drug molecules up to large biopharmaceuticals.



Find in our brochure **Diagnostic Tools** plenty of substrates for reporter enzymes and drug interaction studies, metabolites, glucuronides and inhibitors, inducers, antibody conjugates and cross-linkers, natural products with biological and pharmacological activity, carbohydrates, dyes and fluorescent labels as Tools in Immunology, Diagnostic, Biochemistry and Molecular Biology.



DRUG DISCOVERY - DRUG DELIVERY - DIAGNOSTICS

**Reagents & Resins for Solid Phase Chemistry
Natural & Unusual Amino Acids & Building Blocks**

**The Worldwide largest Selection of Reagents for Drug Delivery
Poly(Amino Acids) - PEGylation - Dendrimers**

**Natural Products, Enzymes & Enzyme Substrates,
Carbohydrates & Fluorophores**

for

Molecular Biology, Microbiology, Biochemistry & Diagnostics

Custom Synthesis of Pharmaceutical Intermediates & Building Blocks

**IRIS BIOTECH GMBH
ADALBERT-ZOELLNER-STR. 1
D-95615 MARKTREDWITZ, GERMANY
PHONE +49 9231 97121 - 0
FAX +49 9231 97121 - 99
INFO@IRIS-BIOTECH.DE
WWW.IRIS-BIOTECH.DE**