


Multifunctional X-ray Film Reader

Digitize Your Traditional Dental X-ray Film



Functions

- ◆ Combines 2.5 inch LCD and X-ray Reader together, image displays on 2.5 inch LCD directly.
- ◆ Can connect to our intra oral camera either by wired or wireless.
- ◆ Multifunctional, can switch X-ray and intra oral camera images anytime.
- ◆ Built-in rechargeable battery. Up to 3 hours use when fully charged.
- ◆ Built-in SD card stores up to 500 pictures.
- ◆ Light weight, easy to operate.
- ◆ VIDEO output, high resolution, high quality.
- ◆ PAL and NTSC.



Built-in SD Card



Connect PC



Connect Wired Intraoral Camera



Connect Wireless Intraoral Camera

Specifications

- ◆ Voltage: 8.5V
- ◆ Power: 500 mA
- ◆ Signal Output: Video
- ◆ Dimensions: 103 mm * 75 mm * 95 mm
- ◆ Weight: 550 g
- ◆ Accessories: X-ray salver, Video cable, charger