

# X-ray Film Reader

Digitize Your Traditional Dental X-ray Film


## Features

- ◆ Plug and Play
- ◆ USB connection
- ◆ Pocket Size
- ◆ Dental Software
- ◆ Clear Image
- ◆ User friendly
- ◆ Better communication
- ◆ Affordable



## Functions

- ✦ Captures traditional periodontitis X-ray films to digital files quickly. Save to computer for recording patients' data.
- ✦ Can zoom X-ray film images up to 50 times using included dental software (Windows™).
- ✦ Dentists gain better communication with patients.
- ✦ Easy Installation: Simply connect to computer with USB cable (included).
- ✦ User friendly design: Simple to use by all the dental practitioners and assistants.
- ✦ Image is captured by latest technology CMOS high resolution camera and bright LED lamps.
- ✦ Adjustable light source enables dental practitioners adjust image brightness.
- ✦ Compact design and light weight offers space-saving and portability.



## Specifications

- ✦ Voltage: DC 5V (powered via USB cable)
- ✦ Power: 300 mA
- ✦ Signal output: USB 2.0
- ✦ Dimensions: 70 mm \* 84 mm \* 80 mm
- ✦ Weight: 430 g
- ✦ Accessories: X-ray salver, image software, USB cable