

Datasheet

Product code	:	DRP1380
Category	:	Diabetes Research Peptides
Name	:	(Glu ⁹)-Exenatide (2-39) trifluoroacetate salt
Sequence	:	H-Gly-Glu-Gly-Thr-Phe-Thr-Ser-Glu-Leu-Ser-Lys-Gln-Met-Glu-Glu-Glu-Ala-Val-Arg-Leu-Phe-Ile-Glu-Trp-Leu-Lys-Asn-Gly-Gly-Pro-Ser-Ser-Gly-Ala-Pro-Pro-Pro-Ser-NH ₂ trifluoroacetate salt
One Letter Code	:	GEGTFTSELSKQMEEEAVRLFIEWLKNGGPSSGAPPPS
Formula	:	C ₁₇₉ H ₂₇₇ N ₄₇ O ₅₉ S
Molecular weight	:	4063.52 gm/mole
CAS	:	[196109-34-9] net
Synonyms	:	(Glu ⁹)-Exendin-4 (2-39), (Des-His ¹ ,Glu ⁹)-Exenatide, (Des-His ¹ ,Glu ⁹)-Exendin-4
Details	:	Exendin-4 analog, effective GLP-1 antagonist.
Specification	:	Purity: min. 95% (RP-HPLC, 214nm) Identity: MALDI-MS Counter ion: trifluoroacetate
Storage	:	-20 ± 5 °C
Price	:	On request