

Datasheet

Product code	:	DRP1440
Category	:	Diabetes Research Peptides
Name	:	(D-Ala ²)-Gastric Inhibitory Polypeptide (human) trifluoroacetate salt
Sequence	:	H-Tyr-D-Ala-Glu-Gly-Thr-Phe-Ile-Ser-Asp-Tyr-Ser-Ile-Ala-Met-Asp-Lys-Ile-His-Gln-Gln-Asp-Phe-Val-Asn-Trp-Leu-Leu-Ala-Gln-Lys-Gly-Lys-Lys-Asn-Asp-Trp-Lys-His-Asn-Ile-Thr-Gln-OH trifluoroacetate salt
One Letter Code	:	YAEGTFISDYSIAMDKIHQQDFVNWLLAQKGKKNDWKHNITQ
formula	:	C ₂₂₆ H ₃₃₈ N ₆₀ O ₆₆ S
molecular weight	:	4983.53 gm/mol
CAS	:	444073-04-5
Synonyms	:	(D-Ala ²)-GIP (human), [DAla2] Gastric Inhibitory Peptide (GIP), human
Details	:	GIP receptor agonist.
Specification	:	Purity: min. 95% (RP-HPLC, 214nm) Identity: MALDI-MS Counter ion: trifluoroacetate
Storage	:	-20 ± 5 °C
Price	:	On request